

Cobolt FIC-03

Mounting plate for laser-to-fiber coupling with beam access




Illustration of possible configuration. Laser and fiber coupling option should be ordered separately.

Cobolt model number	FIC-03
Material	Aluminium
FIC-03 is suitable for use with all Cobolt lasers in the following model:	Cobolt Dual Combiner
dimensions [mm] [inches]	260 x 140 x 51 10.2 x 5.5 x 2.0

Fiber coupling options are available for all Cobolt lasers. The option FIC-03 is recommended for the Cobolt Dual Combiner. This option consists of a mounting plate and a coupling block to mount the fiber coupler onto. Fiber coupling options for either multimode fiber (M/M) or single mode polarization maintaining (SM/PM) fiber can be purchased separately. The FIC-03 is designed to be bolted to a standard metric or imperial optical table.

